

## The most comprehensive and accurate newborn screening test that gives early insight into baby's health

### What is newborn screening?

Almost every child is screened for certain medical conditions shortly after birth. Commonly known as newborn screening, this public health program helps identify babies who may be at an increased risk for serious, but treatable, inherited conditions soon after birth or in childhood.

Most babies who are born with these conditions appear healthy at first. Without early screening, a condition may not be discovered until symptoms appear and it may be too late to prevent serious health consequences. Proactive screening can help parents and pediatricians to identify conditions early in life, so that treatment can be started immediately.

Newborn screening is currently performed only for selected conditions which can be analyzed by biochemical assays from blood, and for which treatment options exist<sup>1</sup>. Each country determines which conditions babies will be screened for, so different countries may test for different conditions. In general, most countries screen for 34 conditions that are recommended. This is also known as the Recommended Uniform Screening Panel (RUSP)<sup>2,3</sup>. Life Baby Test is a more extensive newborn screen. It screens for every condition on the RUSP and all state panels, as well as 188 additional conditions that can benefit from early detection.

The newborn screening covered by the National Health System (NHS) primarily uses biochemical analysis to look at analytes and markers in the blood. Life Baby Test is a genetic test that analyzes DNA to check for variations, or changes that may cause certain diseases. Since many diseases cannot be detected via biochemical blood analysis, Life Baby Test is able to screen for a wider range of conditions than current NHS newborn screening.

### Life Baby Test

Life Baby Test is a highly accurate newborn screening test that analyses baby's DNA for more than **200 conditions**, as well as providing personalized genetic information on the metabolism of more than **30 medications**.

Life Baby Test is a **genetic test** and it only screens for conditions that can be treated with medication, dietary modification, or other therapies.

Early screening with Life Baby Test can help parents and pediatrician to know whether to take proactive steps, sooner, to care for baby's health.

There are two types of Test:

**Standard** – It analyses 86 genes for 87 diseases (including 34 conditions of the Recommended Uniform Screening Panel RUSP) and 12 genes for more than 30 medications.

**Extended** – It analyses 220 genes for 222 diseases (including 34 conditions of the Recommended Uniform Screening Panel RUSP) and 12 genes for more than 30 medications.

## Indications relative to Life Baby Test

- Individuals presenting with the most common symptoms of a metabolic disease, typically include neurological symptoms such as developmental delay, seizure, lethargy, ataxia, behavioural abnormalities, deafness, blindness, and additionally organomegaly and ophthalmologic findings.
- Individuals with a positive family history of metabolic disease<sup>4,5</sup>
- Individuals without a positive family history but with symptoms resembling the specific disease indication
- Individuals with a negative, but suspected, family history in order to perform proper genetic counselling (prenatal analyses are recommended in families of affected individuals)

## What is the procedure for Life Baby Test?

The DNA isolated from cells collected by a buccal swab is then amplified by PCR. Through a state-of-the-art technological process, named massive parallel sequencing (MPS), which uses Next Generation Sequencing (NGS) techniques with ILLUMINA/Thermo Fisher sequencing instruments, genes included in Life Baby panels are completely sequenced (exons and adjacent intronic regions,  $\pm 5$  nucleotides) (Table 1 and 2) at high read depth. The resulting genetic sequences are analysed via an advanced bioinformatics analysis, to check the presence of potential mutations in the genes under investigation.

## Results of Life Baby Test

### “POSITIVE“

#### **Pathogenic / Likely Pathogenic mutation(s) detected**

Indicates that a well characterized disease-causing mutation(s) was identified. This result can indicate the best way to treat the disease.

Only known pathogenic and likely pathogenic mutations are reported

### “NEGATIVE“

#### **NO Pathogenic / Likely Pathogenic mutation(s) detected**

Indicates that no disease causing mutations have been detected in the target genes screened.

A negative test result reduces, but doesn't eliminate, a child's risk of being affected by these genetic conditions or any other genetic conditions not covered by this test.

Life Baby test is only intended to identify mutations in DNA that are very likely to cause genetic conditions (pathogenic or likely pathogenic). Variants of uncertain clinical significance (VOUS), i.e. findings with insufficient evidence available for unequivocal determination of clinical significance, are not reported with Life Baby test.

No single genetic test can detect all of the possible gene variants that could cause a condition. Parental DNA is only analyzed using targeted testing if indicated by the child's test results. A separate report will not be issued to the parent.

### Parameters used to report the genetic variations

The test analyses only the genes listed in Table 1 and 2. Only mutations classified as "known pathogenic" in accordance with the relevant scientific literature and the current classification in the Human Gene Mutation Database (HGMD), updated on the date of the sample collection, will be reported. Moreover, in compliance with the indications of the American College of Medical Genetics (ACMG)<sup>6</sup>, only mutations with a Minor Allele Frequency (MAF) <5% (1000 Genomes Project) are considered as pathogenic or possibly pathogenic; this measurement refers to the frequency in which the less common allele is present in the general population.

### Target Coverage

Life Baby test is a highly accurate newborn screening test that gives early insight into a baby's health. Life Baby test uses an advanced DNA sequencing technology, named Next Generation Sequencing (NGS), coupled with a leading bioinformatic analysis, to detect mutations in 220 genes associated with 222 genetic disorders, with > 99% accuracy. Many of these conditions cannot be detected by carrier screening, standard prenatal tests, or state newborn screening.

### Accuracy of Life Baby Test

Present DNA sequencing techniques are more than 99% accurate. Even though this test is very accurate, the limitations of this examination are always to be taken into consideration. Please read below.

### Limits of Life Baby Test

This examination analyses only genetic diseases and genes listed in Table 1 and 2. The test does not detect other genetic diseases or genes that were not specifically targeted.

Moreover, the test cannot detect:

- Mutation located into intronic regions beyond +/- nucleotides from the breakpoints.
- Deletions, inversions or duplications of more than 20 bps.
- Germline mosaicism (i.e. mutations occurring only in the gametes)

The intrinsic limitation of the NGS methodology is the lack of coverage uniformity of each examined genetic region. Quantity and quality of the DNA extracted from samples is one of the potential causes of such lack of uniformity, which may lead to the lack of detection of gene mutations. Due to this limitation, NGS tests may not detect specific genetic mutations in the selected genes.

Some of these variations may not be identified or validated yet by the scientific community and, therefore, may not be classified as pathogenic variations at the time of analysis. For a correct interpretation of results, we need to have accurate information on the health of the patient and any pathology in the clinical history of the couple and their relatives. This information allows our geneticists to have a better interpretation of genetic results.

## Newborn screening with unmatched detection of serious disorders

The true goal of newborn screening is to detect serious, prevalent, and clinically-actionable diseases. That's why the Life Baby test has been methodically designed by the genetic experts to maximize detection rates for the diseases that matter the most.

The Life Baby test screens for the most clinically relevant and impactful genetic conditions that typically affect health in infancy or childhood. The disorders included in the Life Baby test panel have been carefully selected based on carrier rate, clinical severity, and availability of treatment options. This test only includes conditions that may be treated with medication, dietary modification, or other therapies, so proactive steps can be taken with a pediatrician to care for the child's health.

Every Life Baby test also includes an additional pharmacogenetic analysis of a child's response to more than 30 medications, thus allowing a personalized treatment throughout life. This test is offered at no additional cost and can be completed using the same DNA sample collected for Life Baby.

## List of genes and genetic disorders screened

Life Baby includes genes to cover all ACMG core panel phenotypes for newborn screening

**Table 1:** Genes analysed in the Life Baby **STANDARD** and **EXTENDED** Panel

<b>ABCD1</b>	Adrenoleukodystrophy	<b>Standard Panel</b>
<b>ACAD8</b>	Isobutyryl-CoA dehydrogenase deficiency	<b>Standard Panel</b>
<b>ACADM</b>	Acyl-CoA dehydrogenase, medium chain, deficiency of	<b>Standard Panel</b>
<b>ACADS</b>	Acyl-CoA dehydrogenase, short-chain, deficiency of	<b>Standard Panel</b>
<b>ACADSB</b>	2-methylbutyrylglucosuria	<b>Standard Panel</b>
<b>ACADVL</b>	very long-chain acyl-CoA dehydrogenase deficiency	<b>Standard Panel</b>
<b>ACAT1</b>	Alpha-methylacetoacetic aciduria	<b>Standard Panel</b>
<b>ADA</b>	Severe combined immunodeficiency due to ADA deficiency	<b>Standard Panel</b>
<b>ADK</b>	Hypermethioninemia due to adenosine kinase deficiency	<b>Standard Panel</b>
<b>AHCY</b>	Hypermethioninemia with deficiency of S-adenosylhomocysteine hydrolase	<b>Standard Panel</b>
<b>ARG1</b>	Argininemia	<b>Standard Panel</b>
<b>ASL</b>	Argininosuccinic aciduria	<b>Standard Panel</b>
<b>ASS1</b>	Citrullinemia Type 1	<b>Standard Panel</b>
<b>AUH</b>	3-methylglutaconic aciduria, type I	<b>Standard Panel</b>
<b>BCKDHA</b>	Maple syrup urine disease, type Ia	<b>Standard Panel</b>

<b>BCKDHB</b>	Maple syrup urine disease, type Ib	<b>Standard Panel</b>
<b>BTD</b>	Biotinidase deficiency	<b>Standard Panel</b>
<b>CBS</b>	Homocystinuria, B6-responsive and nonresponsive types	<b>Standard Panel</b>
<b>CFTR</b>	Cystic fibrosis	<b>Standard Panel</b>
<b>CPS1</b>	Carbamoylphosphate synthetase I deficiency	<b>Standard Panel</b>
<b>CPT1A</b>	Carnitine palmitoyltransferase type I deficiency	<b>Standard Panel</b>
<b>CPT2</b>	Carnitine palmitoyltransferase type II deficiency	<b>Standard Panel</b>
<b>CYP21A2</b>	Adrenal hyperplasia, congenital, due to 21-hydroxylase deficiency	<b>Standard Panel</b>
<b>DBT</b>	Maple syrup urine disease, type II	<b>Standard Panel</b>
<b>DNAJC19</b>	Hyperphenylalaninemia, mild, non-BH4-deficient	<b>Standard Panel</b>
<b>ETFA</b>	Glutaric acidemia IIA	<b>Standard Panel</b>
<b>ETFB</b>	Glutaric acidemia IIB	<b>Standard Panel</b>
<b>ETFDH</b>	Glutaric acidemia IIC	<b>Standard Panel</b>
<b>FAH</b>	Tyrosinemia, type I	<b>Standard Panel</b>
<b>G6PD</b>	Hemolytic anemia, G6PD deficient (favism)	<b>Standard Panel</b>
<b>GAA</b>	Glycogen storage disease II - Pompe disease	<b>Standard Panel</b>
<b>GALC</b>	Krabbe disease	<b>Standard Panel</b>
<b>GALE</b>	Galactose epimerase deficiency	<b>Standard Panel</b>
<b>GALK1</b>	Galactokinase deficiency with cataracts	<b>Standard Panel</b>
<b>GALT</b>	Galactosemia	<b>Standard Panel</b>
<b>GBA</b>	Gaucher disease, type I	<b>Standard Panel</b>
<b>GCDH</b>	Glutaricaciduria, type I	<b>Standard Panel</b>
<b>GCH1</b>	Hyperphenylalaninemia, BH4-deficient, B	<b>Standard Panel</b>
<b>GJB2</b>	Deafness, autosomal recessive 1A	<b>Standard Panel</b>
<b>GJB3</b>	Deafness, digenic, GJB2/GJB3	<b>Standard Panel</b>
<b>GJB6</b>	Deafness, digenic GJB2/GJB6	<b>Standard Panel</b>
<b>GLA</b>	Fabry disease	<b>Standard Panel</b>
<b>GSS</b>	Glutathione synthetase deficiency - 5-oxoprolinuria	<b>Standard Panel</b>
<b>HADH</b>	Short-chain hydroxyacyl-coenzyme A dehydrogenase deficiency	<b>Standard Panel</b>

<b>HADHA</b>	long-chain hydroxyacyl-CoA dehydrogenase deficiency	<b>Standard Panel</b>
<b>HADHB</b>	Trifunctional protein deficiency	<b>Standard Panel</b>
<b>HBA1</b>	Thalassemia, alpha-	<b>Standard Panel</b>
<b>HBA2</b>	Thalassemia, alpha-	<b>Standard Panel</b>
<b>HBB</b>	Sickle cell anemia	<b>Standard Panel</b>
<b>HLCS</b>	Holocarboxylase synthetase deficiency	<b>Standard Panel</b>
<b>HMGCL</b>	HMG-CoA lyase deficiency	<b>Standard Panel</b>
<b>HPD</b>	Tyrosinemia, type III	<b>Standard Panel</b>
<b>IDUA</b>	Mucopolysaccharidosis type I <sub>h</sub>	<b>Standard Panel</b>
<b>IL2RG</b>	Severe combined immunodeficiency, X-linked	<b>Standard Panel</b>
<b>IVD</b>	Isovaleric acidemia	<b>Standard Panel</b>
<b>LMBRD1</b>	Methylmalonic aciduria and homocystinuria, cblF type	<b>Standard Panel</b>
<b>MAT1A</b>	Hypermethioninemia, due to methionine adenosyltransferase I/III deficiency	<b>Standard Panel</b>
<b>MCCC1</b>	3-Methylcrotonyl-CoA carboxylase 1 deficiency	<b>Standard Panel</b>
<b>MCCC2</b>	3-Methylcrotonyl-CoA carboxylase 2 deficiency	<b>Standard Panel</b>
<b>MLYCD</b>	Malonyl-CoA decarboxylase deficiency	<b>Standard Panel</b>
<b>MMAA</b>	Methylmalonic aciduria, vitamin B12-responsive	<b>Standard Panel</b>
<b>MMAB</b>	Methylmalonic aciduria, vitamin B12-responsive, due to defect in synthesis of adenosylcobalamin, cblB complementation type	<b>Standard Panel</b>
<b>MMACHC</b>	Methylmalonic aciduria and homocystinuria, cblC type	<b>Standard Panel</b>
<b>MMADHC</b>	Methylmalonic aciduria and homocystinuria, cblD type	<b>Standard Panel</b>
<b>MUT</b>	Methylmalonic aciduria, mut(0) type	<b>Standard Panel</b>
<b>NADK2</b>	2,4-dienoyl-CoA reductase deficiency	<b>Standard Panel</b>
<b>NAGS</b>	N-acetylglutamate synthase deficiency	<b>Standard Panel</b>
<b>OAT</b>	Gyrate atrophy of choroid and retina with or without ornithinemia	<b>Standard Panel</b>
<b>OPA3</b>	3-methylglutaconic aciduria, type III	<b>Standard Panel</b>
<b>OTC</b>	Ornithine transcarbamylase deficiency	<b>Standard Panel</b>
<b>PAH</b>	Phenylketonuria	<b>Standard Panel</b>
<b>PAX8</b>	Hypothyroidism, congenital, due to thyroid dysgenesis or hypoplasia	<b>Standard Panel</b>

<b>PCBD1</b>	Hyperphenylalaninemia, BH4-deficient, D	<b>Standard Panel</b>
<b>PCCA</b>	Propionic acidemia	<b>Standard Panel</b>
<b>PCCB</b>	Propionic acidemia	<b>Standard Panel</b>
<b>PTS</b>	Hyperphenylalaninemia, BH4-deficient, A	<b>Standard Panel</b>
<b>QDPR</b>	Hyperphenylalaninemia, BH4-deficient, C	<b>Standard Panel</b>
<b>SLC22A5</b>	Carnitine deficiency, systemic primary	<b>Standard Panel</b>
<b>SLC25A13</b>	Citrullinemia, type II, adult-onset - neonatal-onset	<b>Standard Panel</b>
<b>SLC25A15</b>	Hyperornithinemia-hyperammonemia-homocitrullinemia syndrome	<b>Standard Panel</b>
<b>SLC25A20</b>	Carnitine-acylcarnitine translocase deficiency	<b>Standard Panel</b>
<b>SLC37A4</b>	Glycogen storage disease Ib	<b>Standard Panel</b>
<b>SMPD1</b>	Niemann-Pick disease, type A	<b>Standard Panel</b>
<b>TAT</b>	Tyrosinemia, type II	<b>Standard Panel</b>
<b>TAZ</b>	3-methylglutaconic aciduria, type II - Barth syndrome	<b>Standard Panel</b>
<b>TSHR</b>	Hypothyroidism, congenital, nongoitrous, 1	<b>Standard Panel</b>
<b>ABCC8</b>	Familial hyperinsulinism ABCC8-related	<b>Extended panel</b>
<b>ABCD1</b>	Adrenoleukodystrophy	<b>Extended panel</b>
<b>ABCD4</b>	Methylmalonic aciduria and homocystinuria, cblJ type	<b>Extended panel</b>
<b>ACAD8</b>	Isobutyryl-CoA dehydrogenase deficiency	<b>Extended panel</b>
<b>ACAD9</b>	acyl-CoA dehydrogenase-9 (ACAD9) deficiency	<b>Extended panel</b>
<b>ACADM</b>	Acyl-CoA dehydrogenase, medium chain, deficiency of	<b>Extended panel</b>
<b>ACADS</b>	Acyl-CoA dehydrogenase, short-chain, deficiency of	<b>Extended panel</b>
<b>ACADSB</b>	2-methylbutyrylglucosaminuria	<b>Extended panel</b>
<b>ACADVL</b>	very long-chain acyl-CoA dehydrogenase deficiency	<b>Extended panel</b>
<b>ACAT1</b>	Alpha-methylacetoacetic aciduria	<b>Extended panel</b>
<b>ACSF3</b>	Combined malonic and methylmalonic aciduria	<b>Extended panel</b>
<b>ADA</b>	Severe combined immunodeficiency due to ADA deficiency	<b>Extended panel</b>
<b>ADK</b>	Hypermethioninemia due to adenosine kinase deficiency	<b>Extended panel</b>
<b>AGL</b>	Glycogen storage disease, type III	<b>Extended panel</b>
<b>AGXT</b>	Primary hyperoxaluria, type 1	<b>Extended panel</b>

<b>AHCY</b>	Hypermethioninemia with deficiency of S-adenosylhomocysteine hydrolase	<b>Extended panel</b>
<b>AK2</b>	Reticular dysgenesis	<b>Extended panel</b>
<b>AKR1D1</b>	Bile acid synthesis defect, congenital, 3	<b>Extended panel</b>
<b>ALDH4A1</b>	Hyperprolinemia, type II	<b>Extended panel</b>
<b>ALDH7A1</b>	Epilepsy, pyridoxine-dependent	<b>Extended panel</b>
<b>ALDOB</b>	Fructose intolerance, hereditary	<b>Extended panel</b>
<b>ALPL</b>	Hypophosphatasia	<b>Extended panel</b>
<b>ANK1</b>	Spherocytosis, type 2	<b>Extended panel</b>
<b>AQP2</b>	Diabetes insipidus, nephrogenic	<b>Extended panel</b>
<b>ARG1</b>	Argininemia	<b>Extended panel</b>
<b>ARSA</b>	Metachromatic leukodystrophy	<b>Extended panel</b>
<b>ARSB</b>	Mucopolysaccharidosis type VI (Maroteaux-Lamy)	<b>Extended panel</b>
<b>ASL</b>	Argininosuccinic aciduria	<b>Extended panel</b>
<b>ASS1</b>	Citrullinemia Type 1	<b>Extended panel</b>
<b>AUH</b>	3-methylglutaconic aciduria, type I	<b>Extended panel</b>
<b>AVPR2</b>	Nephrogenic syndrome of inappropriate antidiuresis / Nephrogenic diabetes insipidus AVPR2-related	<b>Extended panel</b>
<b>BCKDHA</b>	Maple syrup urine disease, type Ia	<b>Extended panel</b>
<b>BCKDHB</b>	Maple syrup urine disease, type Ib	<b>Extended panel</b>
<b>BTD</b>	Biotinidase deficiency	<b>Extended panel</b>
<b>BTK</b>	Agammaglobulinemia, X-linked 1	<b>Extended panel</b>
<b>CASR</b>	Neonatal hyperparathyroidism / Autosomal dominant hypocalcemia	<b>Extended panel</b>
<b>CBS</b>	Homocystinuria, B6-responsive and nonresponsive types	<b>Extended panel</b>
<b>CD247</b>	Immunodeficiency 25	<b>Extended panel</b>
<b>CD320</b>	Methylmalonic aciduria, transient, due to transcobalamin receptor defect	<b>Extended panel</b>
<b>CD3D</b>	Immunodeficiency 19	<b>Extended panel</b>
<b>CD3E</b>	Immunodeficiency 18	<b>Extended panel</b>
<b>CFTR</b>	Cystic fibrosis	<b>Extended panel</b>
<b>COL4A3</b>	Alport syndrome COL4A3-related	<b>Extended panel</b>

<b>COL4A4</b>	Alport syndrome, autosomal recessive	<b>Extended panel</b>
<b>COL4A5</b>	Alport syndrome	<b>Extended panel</b>
<b>CPS1</b>	Carbamoylphosphate synthetase I deficiency	<b>Extended panel</b>
<b>CPT1A</b>	Carnitine palmitoyltransferase type I deficiency	<b>Extended panel</b>
<b>CPT2</b>	Carnitine palmitoyltransferase type II deficiency	<b>Extended panel</b>
<b>CTH</b>	Cystathioninuria	<b>Extended panel</b>
<b>CTNS</b>	Cystinosis	<b>Extended panel</b>
<b>CYBA</b>	Chronic granulomatous disease, autosomal, due to deficiency of CYBA	<b>Extended panel</b>
<b>CYBB</b>	Chronic granulomatous disease CYBB-related	<b>Extended panel</b>
<b>CYP11B1</b>	Congenital adrenal hyperplasia due to 11-beta-hydroxylase deficiency	<b>Extended panel</b>
<b>CYP11B2</b>	Corticosterone methyloxidase deficiency	<b>Extended panel</b>
<b>CYP21A2</b>	Adrenal hyperplasia, congenital, due to 21-hydroxylase deficiency	<b>Extended panel</b>
<b>CYP27A1</b>	Cerebrotendinous xanthomatosis	<b>Extended panel</b>
<b>DBT</b>	Maple syrup urine disease, type II	<b>Extended panel</b>
<b>DCLRE1C</b>	Omenn syndrome / Severe combined immunodeficiency, Athabaskan-type	<b>Extended panel</b>
<b>DEC1</b>	2,4-dienoyl-CoA reductase deficiency	<b>Extended panel</b>
<b>DLD</b>	Dihydrolipoamide dehydrogenase deficiency	<b>Extended panel</b>
<b>DNAJC19</b>	Hyperphenylalaninemia, mild, non-BH4-deficient	<b>Extended panel</b>
<b>DUOX2</b>	Thyroid dysmorphogenesis 6	<b>Extended panel</b>
<b>DUOXA2</b>	Thyroid dysmorphogenesis 5	<b>Extended panel</b>
<b>EPB42</b>	Spherocytosis, type 6	<b>Extended panel</b>
<b>ETFA</b>	Glutaric acidemia IIA	<b>Extended panel</b>
<b>ETFB</b>	Glutaric acidemia IIB	<b>Extended panel</b>
<b>ETFDH</b>	Glutaric acidemia IIC	<b>Extended panel</b>
<b>ETHE1</b>	Ethylmalonic encephalopathy	<b>Extended panel</b>
<b>F9</b>	Factor IX deficiency	<b>Extended panel</b>
<b>FAH</b>	Tyrosinemia, type I	<b>Extended panel</b>
<b>FBN1</b>	Marfan syndrome and other FBN1-related disorders	<b>Extended panel</b>
<b>FBP1</b>	Fructose-1,6-bisphosphatase deficiency	<b>Extended panel</b>

<b>FOLR1</b>	Neurodegeneration due to cerebral folate transport deficiency	<b>Extended panel</b>
<b>FTCD</b>	Glutamate formiminotransferase deficiency	<b>Extended panel</b>
<b>G6PC</b>	Glycogen storage disease, type Ia	<b>Extended panel</b>
<b>G6PD</b>	Hemolytic anemia, G6PD deficient (favism)	<b>Extended panel</b>
<b>GAA</b>	Glycogen storage disease II - Pompe disease	<b>Extended panel</b>
<b>GALC</b>	Krabbe disease	<b>Extended panel</b>
<b>GALE</b>	Galactose epimerase deficiency	<b>Extended panel</b>
<b>GALK1</b>	Galactokinase deficiency with cataracts	<b>Extended panel</b>
<b>GALNS</b>	Mucopolysaccharidosis IVA	<b>Extended panel</b>
<b>GALT</b>	Galactosemia	<b>Extended panel</b>
<b>GAMT</b>	Cerebral creatine deficiency syndrome 2	<b>Extended panel</b>
<b>GATM</b>	Cerebral creatine deficiency syndrome	<b>Extended panel</b>
<b>GBA</b>	Gaucher disease, type I	<b>Extended panel</b>
<b>GCDH</b>	Glutaricaciduria, type I	<b>Extended panel</b>
<b>GCH1</b>	Hyperphenylalaninemia, BH4-deficient, B	<b>Extended panel</b>
<b>GJB2</b>	Deafness, autosomal recessive 1A	<b>Extended panel</b>
<b>GJB3</b>	Deafness, digenic, GJB2/GJB3	<b>Extended panel</b>
<b>GJB6</b>	Deafness, digenic GJB2/GJB6	<b>Extended panel</b>
<b>GLA</b>	Fabry disease	<b>Extended panel</b>
<b>GLIS3</b>	Diabetes mellitus, neonatal, with congenital hypothyroidism	<b>Extended panel</b>
<b>GLUD1</b>	congenital hyperinsulinic hyperammonemia (HI/HA) syndrome	<b>Extended panel</b>
<b>GNAS</b>	Pseudohypoparathyroidism Ia	<b>Extended panel</b>
<b>GNMT</b>	Glycine N-methyltransferase deficiency	<b>Extended panel</b>
<b>GRHPR</b>	Hyperoxaluria, primary, type II	<b>Extended panel</b>
<b>GSS</b>	Glutathione synthetase deficiency - 5-oxoprolinuria	<b>Extended panel</b>
<b>GYS2</b>	Glycogen storage disease 0, liver	<b>Extended panel</b>
<b>HADH</b>	Short-chain hydroxyacyl-coenzyme A dehydrogenase deficiency	<b>Extended panel</b>
<b>HADHA</b>	long-chain hydroxyacyl-CoA dehydrogenase deficiency	<b>Extended panel</b>
<b>HADHB</b>	Trifunctional protein deficiency	<b>Extended panel</b>

<b>HAL</b>	Histidinemia	<b>Extended panel</b>
<b>HAX1</b>	Neutropenia, severe congenital 3, autosomal recessive	<b>Extended panel</b>
<b>HBA1</b>	Thalassemia, alpha-	<b>Extended panel</b>
<b>HBA2</b>	Thalassemia, alpha-	<b>Extended panel</b>
<b>HBB</b>	Sickle cell anemia	<b>Extended panel</b>
<b>HGD</b>	Alkaptonuria	<b>Extended panel</b>
<b>HLCS</b>	Holocarboxylase synthetase deficiency	<b>Extended panel</b>
<b>HMGCL</b>	HMG-CoA lyase deficiency	<b>Extended panel</b>
<b>HMGCS2</b>	HMG-CoA synthase-2 deficiency	<b>Extended panel</b>
<b>HOGA1</b>	Hyperoxaluria, primary, type III	<b>Extended panel</b>
<b>HPD</b>	Tyrosinemia, type III	<b>Extended panel</b>
<b>HSD17B10</b>	17-beta-hydroxysteroid dehydrogenase X (HSD10) deficiency	<b>Extended panel</b>
<b>HSD3B2</b>	Adrenal hyperplasia, congenital, due to 3-beta-hydroxysteroid dehydrogenase 2 deficiency	<b>Extended panel</b>
<b>HSD3B7</b>	Bile acid synthesis defect, congenital, 2	<b>Extended panel</b>
<b>IDS</b>	Mucopolysaccharidosis II	<b>Extended panel</b>
<b>IDUA</b>	Mucopolysaccharidosis type I <sub>h</sub>	<b>Extended panel</b>
<b>IL2RG</b>	Severe combined immunodeficiency, X-linked	<b>Extended panel</b>
<b>IL7R</b>	Severe combined immunodeficiency, T-cell negative, B-cell/natural killer cell-positive type	<b>Extended panel</b>
<b>INS</b>	Diabetes mellitus, permanent neonatal	<b>Extended panel</b>
<b>IVD</b>	Isovaleric acidemia	<b>Extended panel</b>
<b>IYD</b>	Thyroid dysmorphogenesis 5	<b>Extended panel</b>
<b>JAG1</b>	Alagille syndrome 1 / Tetralogy of Fallot	<b>Extended panel</b>
<b>JAK3</b>	Severe combined immunodeficiency	<b>Extended panel</b>
<b>KCNJ11</b>	Familial hyperinsulinism	<b>Extended panel</b>
<b>KCNQ2</b>	Early Infantile epileptic encephalopathy 7 / Benign neonatal seizures 2	<b>Extended panel</b>
<b>LDLR</b>	Familial hypercholesterolemia	<b>Extended panel</b>
<b>LHX3</b>	Combined pituitary hormone deficiency 4	<b>Extended panel</b>

<b>LIG4</b>	LIG4 syndrome	<b>Extended panel</b>
<b>LIPA</b>	Wolman disease / Cholesteryl ester storage disease	<b>Extended panel</b>
<b>LMBRD1</b>	Methylmalonic aciduria and homocystinuria, cblF type	<b>Extended panel</b>
<b>LPL</b>	Lipoprotein lipase deficiency	<b>Extended panel</b>
<b>MAT1A</b>	Hypermethioninemia, due to methionine adenosyltransferase I/III deficiency	<b>Extended panel</b>
<b>MCCC1</b>	3-Methylcrotonyl-CoA carboxylase 1 deficiency	<b>Extended panel</b>
<b>MCCC2</b>	3-Methylcrotonyl-CoA carboxylase 2 deficiency	<b>Extended panel</b>
<b>MCEE</b>	Methylmalonyl-CoA epimerase deficiency	<b>Extended panel</b>
<b>MLYCD</b>	Malonyl-CoA decarboxylase deficiency	<b>Extended panel</b>
<b>MMAA</b>	Methylmalonic aciduria, vitamin B12-responsive	<b>Extended panel</b>
<b>MMAB</b>	Methylmalonic aciduria, vitamin B12-responsive, due to defect in synthesis of adenosylcobalamin, cblB complementation type	<b>Extended panel</b>
<b>MMACHC</b>	Methylmalonic aciduria and homocystinuria, cblC type	<b>Extended panel</b>
<b>MMADHC</b>	Methylmalonic aciduria and homocystinuria, cblD type	<b>Extended panel</b>
<b>MPI</b>	Congenital disorder of glycosylation, type Ib	<b>Extended panel</b>
<b>MPL</b>	Congenital amegakaryocytic thrombocytopenia	<b>Extended panel</b>
<b>MTHFR</b>	Homocystinuria due to MTHFR deficiency	<b>Extended panel</b>
<b>MTR</b>	Homocystinuria-megaloblastic anemia, cobalamin G type	<b>Extended panel</b>
<b>MTRR</b>	Homocystinuria, cobalamin E type	<b>Extended panel</b>
<b>MTTP</b>	Abetalipoproteinemia	<b>Extended panel</b>
<b>MUT</b>	Methylmalonic aciduria, mut(0) type	<b>Extended panel</b>
<b>MVK</b>	Mevalonic aciduria	<b>Extended panel</b>
<b>NADK2</b>	2,4-dienoyl-CoA reductase deficiency	<b>Extended panel</b>
<b>NAGS</b>	N-acetylglutamate synthase deficiency	<b>Extended panel</b>
<b>NHEJ1</b>	Severe combined immunodeficiency with microcephaly, growth retardation, and sensitivity to ionizing radiation	<b>Extended panel</b>
<b>NKX2-1</b>	Choreoathetosis, hypothyroidism, and neonatal respiratory distress	<b>Extended panel</b>
<b>NKX2-5</b>	Hypothyroidism, congenital nongoitrous, 5	<b>Extended panel</b>
<b>NPC1</b>	Niemann-Pick disease, type C1	<b>Extended panel</b>

<b>NPC2</b>	Niemann-pick disease, type C2	<b>Extended panel</b>
<b>OAT</b>	Gyrate atrophy of choroid and retina with or without ornithinemia	<b>Extended panel</b>
<b>OPA3</b>	3-methylglutaconic aciduria, type III	<b>Extended panel</b>
<b>OTC</b>	Ornithine transcarbamylase deficiency	<b>Extended panel</b>
<b>PAH</b>	Phenylketonuria	<b>Extended panel</b>
<b>PAX8</b>	Hypothyroidism, congenital, due to thyroid dysgenesis or hypoplasia	<b>Extended panel</b>
<b>PC</b>	Pyruvate carboxylase deficiency	<b>Extended panel</b>
<b>PCBD1</b>	Hyperphenylalaninemia, BH4-deficient, D	<b>Extended panel</b>
<b>PCCA</b>	Propionic acidemia	<b>Extended panel</b>
<b>PCCB</b>	Propionic acidemia	<b>Extended panel</b>
<b>PHGDH</b>	3-phosphoglycerate dehydrogenase deficiency	<b>Extended panel</b>
<b>PHKB</b>	Glycogen storage disease, type IXb	<b>Extended panel</b>
<b>PNP</b>	Immunodeficiency due to purine nucleoside phosphorylase deficiency	<b>Extended panel</b>
<b>PNPO</b>	Pyridoxamine 5'-phosphate oxidase deficiency	<b>Extended panel</b>
<b>POU1F1</b>	Combined pituitary hormone deficiency 1	<b>Extended panel</b>
<b>PPARG</b>	peroxisome proliferator-activated receptor gamma (PPAR-g) ligand resistance syndrome (PLRS) or familial partial lipodystrophy type 3	<b>Extended panel</b>
<b>PRF1</b>	Hemophagocytic lymphohistiocytosis, familial, 2	<b>Extended panel</b>
<b>PRODH</b>	Hyperprolinemia, type I	<b>Extended panel</b>
<b>PROP1</b>	Combined pituitary hormone deficiency 2	<b>Extended panel</b>
<b>PRRT2</b>	Familial infantile convulsions with paroxysmal choreoathetosis	<b>Extended panel</b>
<b>PTPRC</b>	Severe combined immunodeficiency PTPRC-related	<b>Extended panel</b>
<b>PTS</b>	Hyperphenylalaninemia, BH4-deficient, A	<b>Extended panel</b>
<b>PYGL</b>	Glycogen storage disease, type VI	<b>Extended panel</b>
<b>QDPR</b>	Hyperphenylalaninemia, BH4-deficient, C	<b>Extended panel</b>
<b>RAG1</b>	Omenn syndrome and other RAG1-related disorders	<b>Extended panel</b>
<b>RAG2</b>	Omenn syndrome RAG2-related	<b>Extended panel</b>
<b>RB1</b>	Retinoblastoma	<b>Extended panel</b>
<b>SCN2A</b>	Early infantile epileptic encephalopathy 11 / Benign familial infantile seizures 3	<b>Extended panel</b>

<b>SCN8A</b>	Early infantile epileptic encephalopathy 13 / Benign familial infantile seizures 5	<b>Extended panel</b>
<b>SLC22A5</b>	Carnitine deficiency, systemic primary	<b>Extended panel</b>
<b>SLC25A13</b>	Citrullinemia, type II, adult-onset - neonatal-onset	<b>Extended panel</b>
<b>SLC25A15</b>	Hyperornithinemia-hyperammonemia-homocitrullinemia syndrome	<b>Extended panel</b>
<b>SLC25A20</b>	Carnitine-acylcarnitine translocase deficiency	<b>Extended panel</b>
<b>SLC26A4</b>	Pendred syndrome	<b>Extended panel</b>
<b>SLC2A1</b>	Glucose transporter 1 deficiency syndrome and other SLC2A1-related disorders	<b>Extended panel</b>
<b>SLC3A1</b>	Cystinuria	<b>Extended panel</b>
<b>SLC37A4</b>	Glycogen storage disease Ib	<b>Extended panel</b>
<b>SLC39A4</b>	Acrodermatitis enteropathica	<b>Extended panel</b>
<b>SLC4A1</b>	Distal renal tubular acidosis and other SLC4A1-related disorders	<b>Extended panel</b>
<b>SLC5A5</b>	Thyroid dysmorphogenesis 1	<b>Extended panel</b>
<b>SLC7A7</b>	Lysinuric protein intolerance	<b>Extended panel</b>
<b>SLC7A9</b>	Cystinuria	<b>Extended panel</b>
<b>SMPD1</b>	Niemann-Pick disease, type A	<b>Extended panel</b>
<b>SPR</b>	Sepiapterin reductase deficiency	<b>Extended panel</b>
<b>STAR</b>	Lipoid adrenal hyperplasia	<b>Extended panel</b>
<b>STX11</b>	Hemophagocytic lymphohistiocytosis, familial, 4	<b>Extended panel</b>
<b>SUCLA2</b>	Mitochondrial DNA depletion syndrome 5 (encephalomyopathic with or without methylmalonic aciduria)	<b>Extended panel</b>
<b>SUCLG1</b>	Mitochondrial DNA depletion syndrome 9 (encephalomyopathic type with methylmalonic aciduria)	<b>Extended panel</b>
<b>TAT</b>	Tyrosinemia, type II	<b>Extended panel</b>
<b>TAZ</b>	3-methylglutaconic aciduria, type II - Barth syndrome	<b>Extended panel</b>
<b>TCIRG1</b>	Osteopetrosis 1	<b>Extended panel</b>
<b>TG</b>	Thyroid dysmorphogenesis 4	<b>Extended panel</b>
<b>TH</b>	Segawa syndrome	<b>Extended panel</b>
<b>THRA</b>	Congenital nongoitrous hypothyroidism 6	<b>Extended panel</b>
<b>TMEM70</b>	Mitochondrial complex V (ATP synthase) deficiency, nuclear type 2	<b>Extended panel</b>

<b>TPO</b>	Thyroid dysmorphogenesis 2°	<b>Extended panel</b>
<b>TRHR</b>	Generalized thyrotropin-releasing hormone resistance	<b>Extended panel</b>
<b>TRMU</b>	Acute infantile liver failure	<b>Extended panel</b>
<b>TSHB</b>	Congenital nongoitrous hypothyroidism 4	<b>Extended panel</b>
<b>TSHR</b>	Hypothyroidism, congenital, nongoitrous, 1	<b>Extended panel</b>
<b>TTPA</b>	Ataxia with isolated vitamin E deficiency	<b>Extended panel</b>
<b>UGT1A1</b>	Crigler-Najjar syndrome, types 1 and 2 / Gilbert syndrome	<b>Extended panel</b>
<b>UNC13D</b>	Hemophagocytic lymphohistiocytosis, familial, 3	<b>Extended panel</b>
<b>ZAP70</b>	Immunodeficiency 48	<b>Extended panel</b>
<b>WT1</b>	Wilms tumor, type 1 and other WT1-related disorders	<b>Extended panel</b>

### Additional pharmacogenetic analysis

Every **Life Baby** test includes an additional pharmacogenetic analysis of a child's response to more than **30 medications** that may be prescribed during childhood. This analysis can help doctors to personalize treatments throughout the child's life.

This pharmacogenetic analysis is offered at no additional cost and can be completed using the same DNA sample collected for **Life Baby**.

This test is included with every **Life Baby** test, a genetic screening of **222** childhood genetic conditions.

**Table 2:** List of medications and genes screened

#### **Analgesics**

- Celecoxib (*CYP2C9*)
- Codeine (*CYP2D6*)
- Hydrocodone (*CYP2D6*)
- Oxycodone (*CYP2D6*)
- Tramadol (*CYP2D6*)

#### **Anticoagulants**

- Warfarin (*CYP2C9, VKORC1*)

#### **Anticonvulsant**

- Fosphenytoin (*CYP2C9*)
- Phenytoin (*CYP2C9*)

#### **Antiemetics**

- Ondansetron (*CYP2D6*)

#### **Antifungal Agents**

- Voriconazole (*CYP2C19*)

#### **Antidepressant**

- Amitriptyline (*CYP2D6, CYP2C19*)
- Citalopram (*CYP2C19*)
- Clomipramine (*CYP2D6, CYP2C19*)
- Desipramine (*CYP2D6*)
- Doxepin (*CYP2D6, CYP2C19*)
- Escitalopram (*CYP2C19*)
- Fluoxetine (*CYP2D6*)
- Fluvoxamine (*CYP2D6*)
- Imipramine (*CYP2D6*)
- Nortriptyline (*CYP2D6*)
- Paroxetine (*CYP2D6*)
- Sertraline (*CYP2C19*)
- Trimipramine (*CYP2D6, CYP2C19*)

**Antilipemic Agents**

- Simvastatin (*SLCO1B1*)

**Antineoplastic Agents**

- Capecitabine (*DPYD*)
- Fluorouracil (*DPYD*)
- Mercaptopurine (*TPMT*)
- Thioguanine (*TPMT*)

**Antiplatelet Agents**

- Clopidogrel (*CYP2C19*)

**Antipsychotics**

- Aripiprazole (*CYP2D6*)
- Haloperidone (*CYP2D6*)
- Pimozide (*CYP2D6*)

**Antiretrovirals**

- Atazanavir (*UGT1A1*)

**Immunosuppressant**

- Azathioprine (*TPMT*)
- Tacrolimus (*CYP3A5*)

**Enzyme Inhibitors**

- Eliglustat (*CYP2D6*)

**Psychotropics**

- Atomoxetine (*CYP2D6*)

Please note that almost all individuals will test positive for a variant in at least one of the 12 genes on this panel.

**Pharmacogenetic test reporting.**

Enzymes codified by the genes in table 2 are drug-metabolizing enzymes (DMEs). DMEs are responsible for the drugs metabolic biotransformation and their body elimination. Different relevant polymorphisms have been identified in the gene coding sequences of DMEs enzymes causing loss of therapeutic effects or excessive drug response.

Based on the presence of genetic variants on genes codifying for DMEs enzymes patients can be divided in four phenotypes:

- **Poor Metabolizer - PMs:** PMs have gene inactivating variations on both alleles that result in inactive enzyme or no enzyme at all. PMs have a reduced response or no response and may have increased side effects if treated with standard drug dose.
- **Intermediate Metabolizer – IMs:** IMs have an allele with normal function and a non-functional allele. IMs may require a lower drug dose for a therapeutic effect.
- **Extensive Metabolizer - EMs:** EMs usually have normal drug response as they have two active alleles of the gene.
- **Ultra rapid-Metabolizer - UMs:** UMs usually have an increased expression of genes involved in drug metabolism. UMs may have a reduced or no drug response and may require a higher drug dosage to obtain the therapeutic response. UMs may have three or more active alleles caused by a duplication of the normal allele, or a mutation on one or both alleles causing increased activity.

Eurofins Genoma Group will report the subgroup for the specific drug-metabolizing enzymes identified as polymorphic.

Based on the polymorphisms identified by the pharmacogenetic Life Baby test, the physician can evaluate the adequate drug dosage for the patient. To this purpose, the physician should be aware about metabolic pathway of the specific drug to compare it with the metabolic genetic profile of the patient.

**Polymorphisms investigated:**

*1		None (wild type)	None	Normal function
*2	6	2850 C-T	R296C	Normal function
*3	5	2549delA	frameshift	No function
	3	1749A>G	N166D, frameshift	No function
	5	2549delA		
*4	4	1846G>A	splicing defect	No function
*5		Gene deletion	No protein	No function
*6	3	1707delT	frameshift	No function
*7	6	2935A>C	H324P	No function
*8	3	1758G>T	Stop codon	No function
*9	5	2613-2615delAGA	K281del	Decreased function
*10	1	100C>T	Pro34Ser	Decreased function
*11	intr. 2	883G>C		No function
*12	1	124G>A	Glu42Arg	No function
*14	3	1758G>A	G169R	No function
*17	2	1023C>T	T107I, R296C	Decreased function

	6	2850C>T		
*20	4	1973insG	frameshift	No function
*21	5	2573insC	frameshift	No function
*24	6	2853A>C	I297L	
*38	5	2587-2590delGACT	frameshift	No function
*44	intron 6	2950G>C	splicing defect	No function
*XN		gene amplification	increased protein	Increased function
				Poor metabolizer
*XN (Gene Duplication)				Increased function

**CYP2D6 \*1 and \*2** alleles are associated with a normal enzymatic function, while **CYP2D6 \*3, \*4, \*5, \*6, \*7, \*8, \*9, \*10, \*11, \*12, \*14, \*17, \*20, \*21, \*38, \*44**, correspond to non-functional or decreased activity alleles. Polymorphisms associated to **CYP2D6\*3, CYP2D6\*4, CYP2D6\*5, CYP2D6\*6** alleles are found in 97% of the PM. **CYP2D6\*4** represents the most common allele (21.5-28.6%), followed by the **CYP2D6\*3** (2.7%) and **CYP2D6\*5** (2.6%). Others rare alleles (**CYP2D6\*6, CYP2D6\*7, CYP2D6\*8, CYP2D6\*11, CYP2D6\*12, CYP2D6\*14, CYP2D6\*20, CYP2D6\*31, CYP2D6\*38, CYP2D6\*44**), are associated to non-functional enzyme.

CYP2C19 gene mutations investigated				
CYP2C19 Allele	Exon/intron	Mutation	Effect on protein	Effect on enzymatic function
*1		Nessuna (wild type)	None	Normal function
*2	5	681G>A	splicing defect	No function
*3	4	636G>A	W212X	No function
*4	1	1A>G	start codon mutation	No function
*5	9	1297G>A	R433W	No function
*6	3	395G>A	R132Q	No function
*7	5	1VS5+2T>A	splicing defect	No function
*8	3	358T>C	W120R	No function
*9	3	431G>A	R144H	Decreased function
*10	5	680C>T	P227L	Decreased function
*11	3	449G>A	R150H	
*2, *3, *4, *5, *6, *7, *8, *9, *10, *11				poor metaboliser

**CYP2C19\*1** allele is associated with normal enzymatic function; individuals homozygous for this allele are **extensive metabolizer (EM)**. **CYP2C19 \*2, \*3, \*4, \*5, \*6, \*7, \*8, \*9, \*10, \*11** correspond to non-functional or decreased alleles activity (**poor metabolizer – PM**).

Polymorphisms associated with **CYP2C19\*2, CYP2C19\*3, CYP2C19\*4, CYP2C19\*5, CYP2C19\*6** and **CYP2C19\*7** occurs in 98% of the PM. The most common alleles are **CYP2C19\*2** (75-85% of the Asiatics and 15% of the Europeans and Afro-americans), and **CYP2C19\*3** (6-10% of the Asiatics and rare in Europeans and Afro-americans).

CYP2C9 gene mutations investigated				
CYP2C9 Allele	Exons/intron	Mutation	Effect on protein	Effect on enzymatic function
*1		None (wild type)	None	Normal
*2	3	430C>T	R144C	Decreased function
*3	7	1075A>C	I359L	Decreased function
*4	7	1076T>C	I359T	Decreased function
*5	7	1080C>G	D360E	Decreased function
*6	5	818delA	frameshift	No function
*8	3	449G>A	R150H	Increased function
*9	5	752A>G	H251R	
*10	5	815A>G	E272G	
*11	7	1003C>T	R335W	Decreased function
*14	3	374 G-A	R125H	Decreased function
*18	7	1075A>C	I359L	Decreased function
*2, *3, *4, *5, *6, *11, *14, *18				poor metaboliser

**CYP2C9\*1** allele is associated with normal enzymatic activity; individuals homozygous for this allele are considered to be **extensive metabolizers (EM)**. Alleles **CYP2C9\*2 (R144C)** and **CYP2C9\*3 (I359L)**, were found to be the most common ones in the Caucasian population with frequencies of 8-18% and of 4-10% respectively. The same alleles were found less frequently in Asiatics and Afro-Americans (0.5 – 4%). **CYP2C9\*4** allele was found only in the Japanese population, while **CYP2C9\*5** and **CYP2C9\*6** alleles were found in almost the 2% of the Afro-Americans. Only the 1-2% of European population is homozygous for the CYP2C9\*2 and CYP2C9\*3 alleles.

Individuals carrying **CYP2C9 \*2, \*3, \*4, \*5, \*6, \*11, \*14, \*18** alleles, are found to be poor metabolizers- PM (decreased enzymatic function or no enzymatic function).

Genotype-phenotype correlation for a better patient management:

- Individuals with CYP2C9\*1/\*1 genotype (70% of the Caucasian population) are normal metabolizers.
- Individuals with CYP2C9\*1/ \*2 (16% of the Caucasian population) and \*1/\*3 (10% of the Caucasian population) genotypes are poor metabolizers.
- Individuals with CYP2C9\*2/\*2 (frequency 1%), \*2/\*3 (frequency 1%), \*3/\*3 (frequency 0.3%) genotypes are found to be very poor metabolizers.

**CYP2C9 \*2, CYP2C9 \*3, CYP2C9 \*4, CYP2C9 \*5 and CYP2C9 \*6** alleles are found in almost 98% of PM individuals.

CYP1A2 gene mutations investigated				
CYP1A2 Allele	Exon/intron	Mutation	Effect on protein	Effect on enzymatic function
*1A	1	None (wild type)	None	Normal
*1C	1	g.-3860G>A		Decreased function
*1F	1	g.-164C>A		Elevated inducibility

**CYP1A2\*1A** allele is considered the reference allele to which all variants are compared (normal induction). Two relevant polymorphisms were identified for the **CYP1A2\*1** allele causing **CYP1A2** enzymatic functional alteration: **CYP1A2\*1C** allele, due to **-3860 G>A** point mutation and associated with a decreased metabolic activity; **CYP1A2\*1F** allele, due to **-163 C>A** point mutation, associated with an increased enzymatic induction,

mostly in smokers. CYP1A2 alleles are distributed in the population as follow: **1F/\*1F ~ 46 %**; **\*1A/\*1F ~ 44%**; **\*1A/\*1A ~ 10%**, indicating that the most common phenotype observed is the inducibility increment.

**Some of the medications influenced by the CYP450 system activity**

Generic medication (Brand)	CYP450 genes				Generic medication (Brand)	CYP450 genes			
				2D6		1A2	2C9		2D6
Amitriptyline (Elavil )	✓	✓	✓	✓	diphenhydramine				✓
Aripiprazole (Abilify )				✓	hydroxyzine				✓
Atomoxetine (Strattera )				✓	loratadine (Alavert , Claritin )				✓
Citalopram (Celexa )			✓	✓	terfenadine				✓
Clomipramine (Anafranil )	✓		✓	✓					
Desipramine (Norpramin )				✓	alprenolol				✓
Diazepam (Diazepam Intensol , Valium )			✓	✓	amiodarone	✓	✓	✓	✓
Escitalopram (Lexapro )			✓	✓	atorvastatin		✓		✓
Fluoxetine (Prozac )	✓	✓		✓	brofaramine	✓			
Fluvoxamine (Luvox )	✓			✓	Carvedilol (Coreg )		✓		✓
Haloperidol (Haldol )	✓			✓	cerivastatin			✓	✓
Imipramine (Tofranil )	✓		✓	✓	Flecainide (Tambocor )	✓			✓
Maprotiline (Ludiomil )	✓			✓	Fuvestatin (Lescol )		✓	✓	✓
Nortriptyline (Pamelor , Aventyl )	✓			✓	Irbesartan (Avalide , Avapro )		✓		
Paroxetine (Paxil )	✓			✓	lidocaine	✓			
Perphenazine (Trilafon )	✓			✓	Losartan (Cozaar , Hyzaar )		✓		
Phenytoin (Dilantin )		✓		✓	lovastatin		✓	✓	
Risperidone (Risperdal )				✓	nicardipine				✓
Sertraline (Zoloft )	✓	✓	✓	✓	Metoprolol (Lopressor , Toprol XL )				✓
Thioridazine (Mellaril )	✓			✓	Mexiletine (Mexitil )	✓			✓
Trimipramine (Surmontil )		✓	✓	✓	nicardipine		✓	✓	✓
Venlafaxine (Effexor )				✓	nifedipine			✓	
Zuclopenthixol (Clopixol )				✓	Propafenone (Rythmol )	✓			✓
					propranolol	✓		✓	✓
ciprofloxacin	✓				quinidine		✓		✓
enoxacin	✓				simvastatin		✓	✓	✓
grepafloxacin	✓				Timolol (Betimol , Timoptic , Blocadren , Istalol )				✓
norfloxacin	✓				verapamil	✓			✓
ofloxacin	✓				Warfarin (Coumadin )		✓		
rifampin		✓							
sparfloxacin	✓				acetaminophen	✓			✓
					Celecoxib (Celebrex )		✓		✓
carbamazepine	✓		✓		Codeine				✓
Diazepam (Valium )			✓		diclofenac		✓		✓
felbamate			✓		flurbiprofen		✓		
mephenytoin			✓		hydrocodone				✓
phenobarbitone			✓		ketoprofen		✓		
Phenytoin (Dilantin )		✓	✓		Ibuprofen (Advil , Motrin )		✓		
topiramate			✓		indomethacin		✓		

**Some of the medications influenced by the CYP450 system activity**

Generic medication (Brand)	CYP450 genes				Generic medication (Brand)	CYP450 genes			
<b>Anti-Diabetics</b>					<b>Analgesics/ Anti-inflammatories</b>				
Glimepiride (Amaryl <sup>®</sup> )		✓			meloxicam		✓		
Glipizide (Glucotrol <sup>®</sup> )		✓			methadone		✓	✓	✓
Glyburide (Diabeta <sup>®</sup> , Micronase <sup>®</sup> , Glynase <sup>®</sup> )		✓	✓		Naproxen (Aleve <sup>®</sup> , Naprosyn <sup>®</sup> )	✓	✓		
nateglinide		✓			perhexiline				✓
Pioglitazone (Actos <sup>®</sup> )		✓			phenacetin				✓
Rosiglitazone (Avandia <sup>®</sup> )		✓			piroxicam		✓		
Tolbutamide		✓	✓		propranolol				✓
Chemioterapici					suprofen		✓		
Ondansetron (Zofran <sup>®</sup> )				✓	Oxycodone (Oxycontin <sup>®</sup> )				✓
Tamoxifen (Nolvadex <sup>®</sup> )				✓	Tramadol (Ultram <sup>®</sup> )				✓
Tropisetron				✓	<b>Antacids</b>				
<b>Antihistamines</b>					Esomeprazole (Nexium <sup>®</sup> )			✓	✓
Astemizole				✓	Lansoprazole (Prevacid <sup>®</sup> )			✓	✓
Azelastine				✓	Omeprazole (Prilosec <sup>®</sup> )		✓	✓	✓
chlorpheniramine				✓	Pantoprazole (Protonix <sup>®</sup> )			✓	
<b>HIV Antiviral</b>					Rabepranole (Aciphex <sup>®</sup> )			✓	
					<b>HIV Antiviral</b>				
					amprenavir		✓	✓	✓
					delavirdine	✓	✓	✓	✓
					efavirenz	✓	✓	✓	✓
					indinavir				✓
					nelfinavir			✓	✓
					ritonavir	✓	✓	✓	✓
					saquinavir				✓

CYP3A4 gene mutations investigated				
CYP3A4 Allele	Exon/intron	Mutation	Effect on protein	Effect on enzymatic function
*1A		Wild-type	None	Normal function
*2	7	g.15713T>C	S222P	Decreased function
*6		g.17662-17663insA	frameshift	Premature transcription termination
*8	5	g.13908G>A	R130Q	Decreased function
*11	11	g.21867C>T	T363M	Decreased function
*12	11	g.21896C>T	L373F	Decreased function
*13	11	g.22026C>T	R416L	Decreased function
*16	7	g.15603C>G	T185S	Decreased function
*17	7	g.15615T>C	F189S	Decreased function
*18	10	g.20070T>C	L293P	Increased function
*20		25889_25890insA	488Frameshift	No function
*26		17633C>T	R268X	Decreased function

### Dugs metabolised by CYP3A4

The most important enzymatic substrates, inhibitors and inducers of the CYP3A4 cytochrome are listed in the tables below:

Substrates, inhibitors and inducers of the CYP3A4		
SUBSTRATES	INHIBITORS	INDUCTORS
<b>Antidepressants</b> (Imipramine, amitriptyline, sertraline, venlafaxine, nefazodone) <b>Benzodiazepine</b> (alprazolam, triazolam, midazolam) <b>Antimycotics</b> (ketoconazole, astemizole) <b>Protease inhibitors</b> (ritonavir, indinavir, nelfinavir, saquinavir)	<b>Antidepressant</b> (Nefazodone > fluvoxamine > Fluoxetine > sertraline, paroxetine, venlafaxine)  <b>Azole antifungal</b> (Ketoconazole, itraconazole, fluconazole)	Carbamazepine Dexamethasone Phenobarbital Phenytoin Rifampicin
<b>Others</b> Terfenadine Verapamil Testosterone Theophylline Carbamazepine Cisapride Dexamethasone Erythromycin Ethinyl estradiol Glyburides Cyclosporine Lovastatin	<b>Others</b> Cimetidine Clarithromycin Diltiazem Protease Inhibitors	

### Examples of drug-drug interaction between drugs metabolised by CYP3A4

Interfering drugs	Interaction mechanism	Clinical consequences
Fluoxetine-calcium antagonists	Fluoxetine, by its metabolite norfluoxetine, inhibits the CYP3A4 activity, limiting the calcium-antagonist agents metabolism.	Nausea, flushing, oedema, headache.
Fluoxetine-alprazolam	Fluoxetine, by its metabolite norfluoxetine, inhibits the CYP3A4 activity, which normally metabolise alprazolam. This drug-drug interaction does not affect the activity of triazolam, which is mainly metabolised in the gastrointestinal tract.	Cognitive or psychomotor function impairment.
Astemizole-ketoconazole	Ketoconazole is a strong selective CYP3A4 inhibitor, able to inhibit the astemizole metabolism almost entirely.	Cardio toxicity.

### TPMT gene mutations investigated

TPMT Allele	Exon/intron	Mutation	Effect on protein	Effect on enzymatic function
*1		None (wild type)	None	Normal function
*2		c.288G>C	A80P	Decreased function
*3A		c.460G>A c.719A>G	A154T Y240C	Decreased function
*3B		c.719A>G	Y240C	Decreased function

At least 20 polymorphisms have been identified on the *TPMT* gene, four of them are associated to decrement or loss of enzymatic function. **TPMT\*3A** and **TPMT\*3C** polymorphisms (Poor Metabolizer) are the most common defective alleles (almost 90% of the defective alleles in the European population).

These polymorphisms are transmitted as an autosomal codominant trait and their combination lead to four different alleles associated to different enzymatic activity levels.

### CYP3A5 gene mutations investigated

CYP3A5 Allele	Mutation	Effect on protein	Effect on enzymatic function
CYP3A5*1A	None (wild type)	None	Normal function
CYP3A5*3A	6986A>G; 31611C>T	Splicing defect	Severely decreased function
CYP3A5*3B	3705C>T; 3709_3710insG; 6986A>G; 31611C>T	H30Y; Splicing defect	Severely decreased function
CYP3A5*3C	6986A>G	Splicing defect	Severely decreased function
CYP3A5*3D	6986A>G; 7249T>G	Splicing defect; L82R	Severely decreased function

<b>CYP3A5*3E</b>	6986A>G; 27050A>G; 31611C>T	Splicing defect	<b>Severely decreased function</b>
<b>CYP3A5*3F</b>	6986A>G; 31551T>C; 31611C>T	Splicing defect; I488T	<b>Severely decreased function</b>
<b>CYP3A5*3G</b>	6986A>G; 12952T>C; 31611C>T	Splicing defect	<b>Severely decreased function</b>
<b>CYP3A5*3H</b>	6986A>G; 13108T>C; 31611C>T	Splicing defect	<b>Severely decreased function</b>
<b>CYP3A5*3I</b>	6986A>G; 16903A>G; 31611C>T	Splicing defect	<b>Severely decreased function</b>
<b>CYP3A5*3J</b>	6986A>G; 29782A>G; 31611C>T	Splicing defect; I456V	<b>Severely decreased function</b>
<b>CYP3A5*3K</b>	6986A>G; 29753T>C; 31611C>T	Splicing defect; F446S	<b>Decreased function</b>
<b>CYP3A5*3L</b>	3775A>G; 6986A>G	Y53C; Splicing defect	<b>Decreased function</b>
<b>CYP3A5*6</b>	14690G>A	Splicing defect	<b>No function or severely decreased function</b>
<b>CYP3A5*8</b>	3699C>T	R28C	<b>Decreased function</b>
<b>CYP3A5*9</b>	19386G>A	A337T	<b>Decreased function</b>

CYP3A7 gene mutations investigated				
Gene	Allele	Mutation	Effect on protein	Effect on enzymatic function
CYP3A7	*1C	c.-291G>T -284T>A -282T>C -281A>T -270T>G -262T>A -232A>C		Increased function
	*2	c.1226C>G	T409R	Increased function

UGT1A1 gene mutations investigated				
Allele	Exon/intron	Mutation	Effect on protein	Effect on enzymatic function
*1		promoter repeat [TA] 6	wild type	
*28		promoter repeat [TA] 7		Decreased function
*36		promoter repeat [TA] 5		Increased activity
*37		promoter repeat [TA] 8		Decreased function

Three principal genotypes in the promoter region of the **UGT1A1** gene have been observed in the general population: the wild type homozygous genotype, characterised by 6 thymine-adenine dinucleotide repeats **(TA) 6/6**, the homozygous genotype with **7 TA** dinucleotide repeats **(TA) 7/7**, and the heterozygous **(TA) 6/7** genotype carrying a wild-type allele and the 7 TA expanded allele.

The presence of the 7TA expanded allele lead to a **decreased enzymatic activity**.

## DPYD

Different polymorphisms in the DPYD gene have been found to decrease or abolish the DPD enzymatic function, with the consequent increased risk of severe toxicity, potentially fatal. The most common is the point mutation **G→A** in a splicing site of the exon 14 (**IVS14+ 1G>A**). This nucleotide substitution results in exon loss and the production of a non-functional incomplete protein. Approximately the 1% of the general population carries this allelic variant (known as **DPYD\*2A**) corresponding to almost the 50% of DPD-deficiency. In particular, patients heterozygous for these polymorphisms show partial DPD enzymatic activity, while patients carrying homozygous mutation show DPD enzymatic deficiency. A dose reduction should be considered for patients carrying these variants.

VKORC1 gene mutations investigated		
VKORC1 allele	Mutations	Effect on enzymatic function
-1639	-1639G>A	Decreased Transcription

The **-1639G>A** polymorphism in the promoter of the **VKORC1** (Vitamin K epOxide Reductase Complex subunit 1) gene, target of dicumarolic drugs, is found to be associated to a decreased gene expression. As dicumarolic drugs act by inhibiting Vitamin K epoxide reductase, the presence of the genetic variant (which induce a lower enzyme product resulting in increased drug sensitivity) requires a reduction in drug dose.

## Bibliografia

1. Feuchtbaum, L, et al. Birth prevalence of disorders detectable through newborn screening by race/ethnicity. *Genet. Med.* 2012; 14(11):937-45. PMID: 22766612
2. <https://www.hrsa.gov/advisory-committees/heritable-disorders/rusp/index.html>
3. “Newborn Screening: Beyond the Spot” Urv et Al. *Adv Exp Med Biol.* 2017;1031:323-346.
4. Fuller M, et Al. Epidemiology of lysosomal storage diseases: an overview. In: Mehta A, Beck M, Sunder-Plassmann G, editors. *Fabry Disease: Perspectives from 5 Years of FOS.* Oxford: Oxford PharmaGenesis; 2006.
5. Platt FM, Boland B, van der Spoel AC. Lysosomal storage disorders: the cellular impact of lysosomal dysfunction. *The Journal of cell biology.* 2012 Nov 26;199(5):723-34
6. “ACMG recommendations for reporting of incidental findings in clinical exome and genome sequencing.” Grody et Al. *Genet Med.* 2013 5(7):565-74.